

Urine Analysis Test Strip

UA - 101s & UA -110s

The urine analysis reagent strips are firm plastic strips which several separate reagent area are fixed.

Urine analysis is a useful procedure as an indicator of health or disease, and such, is a part of routine health screenings. It can be used in general evaluation of health, and aids in the diagnosis and monitoring of metabolic or systemic diseases that affect kidney function, endocrine disorders and diseases or disorders of the urinary tract.

Specification :

- Sample : Urine
- Procedure : Easy to use
- Long expiry date
- Package : 50 tests or 100 tests / tube

Urine Analysis Strip

- Ascorbic Acid
- Glucose
- Bilirubin
- Ketones
- Specific Gravity
- Blood
- pH
- Protein
- Urobilinogen
- Nitrite
- Leucocytes

Interpretation of Result

Step 1

Step 2

Step 3

